
## neoligaments"



outcomes for patients, and we want to hear from you.

Help us create the future of sports medicine and share your opinions and expertise with us at: www.neoligaments.com

Distributed in Australia & New Zealand by:




+ 61 (2) 9316 6665 info@neosurgical.com.au

## Why Infinity-Lock Button System?

- Use the Infinity-Lock Button System for AC joint separation requiring a simple and replicable operative solution.
- Device and surgical technique designed by leaders in upper limb surgery.
- One size fits all, eliminating potential for mis-sizing and stocking issues.
- Woven polyester implant utilizes medical manufacturing techniques supported by over 30 years of data.

- Titanium alloy Button provides high strength, low profile fixation in the clavicle.
- Infinity-Lock Button System is a strong and reliable alternative ACJ repair method.


## ordering information



The information above MUST NOT be used to perform the surgical procedure. Always refer to the Instructions for Use LAB 253 which contains detailed information and essential warnings and precautions.

102-1089 Infinity-Lock Button System (supplied sterile) includes:

Infinity-Lock Tube-Tape, 7 mm x 240 mm
Infinity-Lock Button, 4 mm x 12 mm
Cannulated drill bit, plain shank to fit Jacobs chuck, 4.0 mm diameter x 120 mm
Guidewire, diameter 2.0 mm x 150 mm



Optional Disposable Instrument (supplied sterile):

202-1411 CC-Hook, with a curved end, Left 202-1413 CC-Hook, with a curved end, Right



Developed and manufactured by Xiros™
Springfield House, Whitehouse Lane, Leeds LS19 7UE
Tel. +44 (0) 113 238 7202 Fax. +44 (0) 113 238 7201
enquiries@neoligaments.com
www.neoligaments.com

Xiros Limited, registered in England No. 1664824

Neoligaments, Xiros, and Infinity-Lock are trademarks of Xiros. All rights reserved © Neoligaments™ 2019. Worldwide patents and patents pending.